

A.S.A.P.[®]

All Surface Access Polishers



More Shine. Less Time.

Diamond-impregnated **A.S.A.P.[®] Polishers** provide all surface access polishing and are ideal for use on all composite surfaces, lithium disilicate, bisacryls, most ceramics, and zirconia.

- The A.S.A.P. Pre-polisher (purple) is embedded with diamond particles (44 microns) and gently diminishes small surface defects without affecting anatomy, and prepares the surface for a final high gloss shine.
- The A.S.A.P. Final High Shine Polisher (peach) is embedded with small, 3-6 micron sized diamond particles to provide a life-like polish in as little as 20 seconds.
- Autoclavable up to 30 times, A.S.A.P. Polishers can withstand wear while providing an efficient and consistent polish use after use.


All Surface Access Polishing




Proven Performance

A recent university-based study,** comparing the gloss produced by A.S.A.P. Polishers (Clinician's Choice[®]) and Enhance[®]/PoGo[®] (Dentsply Sirona) found that A.S.A.P. Polishers produced a higher gloss results, after 30 seconds, than Enhance and Pogo did after 90 seconds.

Maximum Gloss Per System




Lawson, N. University of Alabama at Birmingham, School of Dentistry. Gloss of composites with the use of two polishing systems. April 2022.
*Enhance[®] and PoGo[®] are registered trademarks of Dentsply Sirona. ** Data on file.



“A.S.A.P. Polishers save chair time by quickly polishing porcelain and composite restorations in **less than a minute.**”

- Danièle Larose, DMD, AAACD





CLINICAL Instant Replay

In a rush?
Watch this quick technique tip!

POWERED BY
Hands-OnLive

Scan to watch

Achieving a Superior Polish with A.S.A.P. All Surface Access Polishers

Bob Margeas, DDS

CLINICAL Instant Replay

In a rush?
Watch this quick technique tip!

POWERED BY
Hands-OnLive

Scan to watch

Composite Finishing Technique for an Esthetic Result

Bob Margeas, DDS



Request a complimentary sample:
sample.clinicianchoice.com/asap-polisher-family



Scan to Read
A.S.A.P. Polishers Dental Advisor Evaluation



| | |
|--|---|
| <p>A.S.A.P. STARTER KIT (130917) Contains: 3 x Pre-polishers Small, 3 x Pre-polishers Large, 3 x Final High Shine Polishers Small, 3 x Final High Shine Polishers Large</p> | <p>A.S.A.P. TRIAL 2-PACK Contains: 1 x Pre-polisher, 1 x Final High Shine Polisher Small (130910) Large (130911)</p> |
| <p>A.S.A.P. REFILL 6-PACK Contains: 3 x Pre-polishers, 3 x Final High Shine Polishers Small (130913) Large (130912)</p> | <p>A.S.A.P. BULK PACK SMALL PRE-POLISHER (130927) Contains: 12 x Pre-polishers Small</p> |
| <p>A.S.A.P. PRE-POLISHER 3-PACK Small (131011) Large (131012)</p> | <p>A.S.A.P. BULK PACK LARGE PRE-POLISHER (130928) Contains: 12 x Pre-polishers Large</p> |
| <p>A.S.A.P. FINAL HIGH SHINE POLISHER 3-PACK Small (132013) Large (132014)</p> | <p>A.S.A.P. BULK PACK SMALL FINAL HIGH SHINE POLISHER (130925) Contains: 12 x Final High Shine Polishers Small</p> |
| | <p>A.S.A.P. BULK PACK LARGE FINAL HIGH SHINE POLISHER (130926) Contains: 12 x Final High Shine Polishers Large</p> |

To view the latest news, education and more, visit:

www.clinicianchoice.com